

Human CRYAA/Alpha-crystallin A chain ELISA Kit

Catalog No: #EK19648

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CRYAA/Alpha-crystallin A chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-crystallin A chain
Target Name	Alpha-crystallin A chain
Other Names	Heat shock protein beta-4,HspB4,CRYA1,HSPB4
Accession No.	P02489
Uniprot	P02489
GeneID	102724652;1409;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.802ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CRYAA

Note: This product is for in vitro research use only