

Mouse Cryaa/Alpha-crystallin A chain ELISA Kit

Catalog No: #EK19649

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cryaa/Alpha-crystallin A chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-crystallin A chain
Target Name	Alpha-crystallin A chain
Other Names	Crya1
Accession No.	P24622
Uniprot	P24622
GeneID	12954;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.45ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cryaa

Note: This product is for in vitro research use only