

Human CPM/Carboxypeptidase M ELISA Kit

Catalog No: #EK19650

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CPM/Carboxypeptidase M ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carboxypeptidase M
Target Name	Carboxypeptidase M
Other Names	CPM
Accession No.	P14384
Uniprot	P14384
GeneID	1368;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CPM

Note: This product is for in vitro research use only