

Mouse Caprin1/Caprin-1 ELISA Kit

Catalog No: #EK19652

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Caprin1/Caprin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caprin-1
Target Name	Caprin-1
Other Names	Cytoplasmic activation- and proliferation-associated protein 1,GPI-anchored membrane protein 1,GPI-anchored protein p137,GPI-p137,Membrane component chromosome 11 surface marker 1,RNA granule protein 105,Gpiap,Gpiap1,Gpip137,M11s1,Rng105,G5E5
Accession No.	Q60865
Uniprot	Q60865
GenID	53872;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Caprin1

Note: This product is for in vitro research use only