

Mouse Btg3/Protein BTG3 ELISA Kit

Catalog No: #EK19653

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Btg3/Protein BTG3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Protein BTG3
Target Name	Protein BTG3
Other Names	Abundant in neuroepithelium area protein,BTG family member 3,Protein Tob5,Ana,Tob5
Accession No.	P50615
Uniprot	P50615
GeneID	12228;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.079ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Btg3

Note: This product is for in vitro research use only