

Human ANXA6/Annexin A6 ELISA Kit

Catalog No: #EK19654

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANXA6/Annexin A6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Annexin A6
Target Name	Annexin A6
Other Names	67 kDa calelectrin, Annexin VI, Annexin-6, Calphobindin-II, CPB-II, Chromobindin-20, Lipocortin VI, Protein III, p68, p70, ANX6
Accession No.	P08133
Uniprot	P08133
GeneID	309;

Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 39.13 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANXA6

Note: This product is for in vitro research use only