

## Mouse Anxa6/Annexin A6 ELISA Kit

Catalog No: #EK19655

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Anxa6/Annexin A6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Annexin A6
Target Name	Annexin A6
Other Names	67 kDa calelectrin, Annexin VI, Annexin-6, Calphobindin-II, CPB-II, Chromobindin-20, Lipocortin VI, Protein III, p68, p70, Anx6
Accession No.	P14824
Uniprot	P14824
GeneID	11749;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 39.14 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Anxa6

Note: This product is for in vitro research use only