

## Rat Anxa6/Annexin A6 ELISA Kit

Catalog No: #EK19656

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Anxa6/Annexin A6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Annexin A6
Target Name	Annexin A6
Other Names	Annexin VI, Annexin-6, Calcium-binding protein 65/67, CBP 65/67, Anx6
Accession No.	P48037
Uniprot	P48037

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.08ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Anxa6

Note: This product is for in vitro research use only