

Rat Actn4/Alpha-actinin-4 ELISA Kit

Catalog No: #EK19658

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Actn4/Alpha-actinin-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-actinin-4
Target Name	Alpha-actinin-4
Other Names	Non-muscle alpha-actinin 4
Accession No.	Q9QXQ0
Uniprot	Q9QXQ0
GeneID	63836;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.46ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Actn4

Note: This product is for in vitro research use only