

Mouse Chi3l1/Chitinase-3-like protein 1 ELISA Kit

Catalog No: #EK19659

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Chi3l1/Chitinase-3-like protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Chitinase-3-like protein 1
Target Name	Chitinase-3-like protein 1
Other Names	BRP39 protein, Breast regression protein 39, Cartilage glycoprotein 39, CGP-39, Brp39, Chil1
Accession No.	Q61362
Uniprot	Q61362
GeneID	12654;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Chi3l1

Note: This product is for in vitro research use only