

Mouse Adam12/Disintegrin and metalloproteinase domain-containing protein 12 ELISA Kit



Catalog No: #EK19661

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adam12/Disintegrin and metalloproteinase domain-containing protein 12 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 12
Target Name	Disintegrin and metalloproteinase domain-containing protein 12
Other Names	ADAM 12,Meltrin-alpha,Mltna
Accession No.	Q61824
Uniprot	Q61824
GeneID	11489;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.47ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adam12

Note: This product is for in vitro research use only