

Mouse Col27a1/Collagen alpha-1(XXVII) chain ELISA Kit



Catalog No: #EK19663

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Col27a1/Collagen alpha-1(XXVII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(XXVII) chain
Target Name	Collagen alpha-1(XXVII) chain
Other Names	Kiaa1870
Accession No.	Q5QNQ9
Uniprot	Q5QNQ9
GeneID	373864;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.081ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Col27a1

Note: This product is for in vitro research use only