

Mouse Col7a1/Collagen alpha-1(VII) chain ELISA Kit

Catalog No: #EK19665

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Col7a1/Collagen alpha-1(VII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collagen alpha-1(VII) chain
Target Name	Collagen alpha-1(VII) chain
Other Names	Long-chain collagen,LC collagen
Accession No.	Q63870
Uniprot	Q63870
GeneID	12836;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.19ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Col7a1

Note: This product is for in vitro research use only