

# Mouse Bmp15/Bone morphogenetic protein 15 ELISA Kit



Catalog No: #EK19666

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bmp15/Bone morphogenetic protein 15 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 15
Target Name	Bone morphogenetic protein 15
Other Names	BMP-15,Growth/differentiation factor 9B,GDF-9B,Gdf9b
Accession No.	Q9Z0L4
Uniprot	Q9Z0L4
GeneID	12155;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.192ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bmp15

Note: This product is for in vitro research use only