

## Mouse Cdh2/Cadherin-2 ELISA Kit

Catalog No: #EK19668

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdh2/Cadherin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-2
Target Name	Cadherin-2
Other Names	Neural cadherin,N-cadherin,CD325
Accession No.	P15116
Uniprot	P15116
GeneID	12558;

## Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.279ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cdh2

Note: This product is for in vitro research use only