

Rat Cdh2/Cadherin-2 ELISA Kit

Catalog No: #EK19669

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cdh2/Cadherin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cadherin-2
Target Name	Cadherin-2
Other Names	Neural cadherin, N-cadherin, CD325
Accession No.	Q9Z1Y3
Uniprot	Q9Z1Y3
GeneID	83501;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ml Sensitivity: 24pg/ml Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdh2

Note: This product is for in vitro research use only