

# Bovine CELA1/Chymotrypsin-like elastase family member 1 ELISA Kit

Catalog No: #EK19671

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CELA1/Chymotrypsin-like elastase family member 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Chymotrypsin-like elastase family member 1
Target Name	Chymotrypsin-like elastase family member 1
Other Names	Elastase I,Elastase-1,ELA1
Accession No.	Q28153
Uniprot	Q28153
GeneID	281139;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.27pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CELA1

Note: This product is for in vitro research use only