

# Pig CELA1/Chymotrypsin-like elastase family member 1 ELISA Kit



Catalog No: #EK19673

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig CELA1/Chymotrypsin-like elastase family member 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Chymotrypsin-like elastase family member 1
Target Name	Chymotrypsin-like elastase family member 1
Other Names	Elastase-1,ELA1
Accession No.	P00772
Uniprot	P00772
GeneID	396766;

## Application Details

Assay Type: Sandwich Detect Range:0.39-25ng/mL Sensitivity:0.121ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CELA1

Note: This product is for in vitro research use only