

Rat Cela1/Chymotrypsin-like elastase family member 1 ELISA Kit



Catalog No: #EK19674

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cela1/Chymotrypsin-like elastase family member 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chymotrypsin-like elastase family member 1
Target Name	Chymotrypsin-like elastase family member 1
Other Names	Elastase-1,Ela1
Accession No.	P00773
Uniprot	P00773
GeneID	24331;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 35pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cela1

Note: This product is for in vitro research use only