

## Rat Calr/Calreticulin ELISA Kit

Catalog No: #EK19676

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Calr/Calreticulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calreticulin
Target Name	Calreticulin
Other Names	CALBP,CRP55,Calcium-binding protein 3,CABP3,Calregulin,Endoplasmic reticulum resident protein 60,ERp60,HACBP
Accession No.	P18418
Uniprot	P18418
GeneID	64202;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.067ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Calr

Note: This product is for in vitro research use only