

# Human CLEC7A/C-type lectin domain family 7 member A ELISA Kit



Catalog No: #EK19677

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CLEC7A/C-type lectin domain family 7 member A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-type lectin domain family 7 member A
Target Name	C-type lectin domain family 7 member A
Other Names	Beta-glucan receptor,C-type lectin superfamily member 12,Dendritic cell-associated C-type lectin 1,DC-associated C-type lectin 1,BGR,CLECSF12,DECTIN1,UNQ539/PRO1082
Accession No.	Q9BXN2
Uniprot	Q9BXN2
GeneID	64581;

## Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:15.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CLEC7A

Note: This product is for in vitro research use only