

Human ADRA2A/Alpha-2A adrenergic receptor ELISA Kit



Catalog No: #EK19678

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ADRA2A/Alpha-2A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-2A adrenergic receptor
Target Name	Alpha-2A adrenergic receptor
Other Names	Alpha-2 adrenergic receptor subtype C10,Alpha-2A adrenoreceptor,Alpha-2A adrenoceptor,ADRA2R,ADRAR
Accession No.	P08913
Uniprot	P08913
GeneID	150;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.079ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ADRA2A

Note: This product is for in vitro research use only