

Rat Adra2a/Alpha-2A adrenergic receptor ELISA Kit

Catalog No: #EK19679

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Adra2a/Alpha-2A adrenergic receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-2A adrenergic receptor
Target Name	Alpha-2A adrenergic receptor
Other Names	Alpha-2A adrenoreceptor, Alpha-2A adrenoceptor, Alpha-2D adrenergic receptor, CA2-47
Accession No.	P22909
Uniprot	P22909
GeneID	25083;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 8.02pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adra2a

Note: This product is for in vitro research use only