

# Mouse Ang2/Angiogenin-2 ELISA Kit

Catalog No: #EK19680

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ang2/Angiogenin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiogenin-2
Target Name	Angiogenin-2
Other Names	Angiogenin-related protein,Angrp
Accession No.	Q64438
Uniprot	Q64438
GeneID	11731;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.17ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ang2

Note: This product is for in vitro research use only