

Mouse C9/Complement component C9 ELISA Kit

Catalog No: #EK19681

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse C9/Complement component C9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement component C9
Target Name	Complement component C9
Accession No.	P06683
Uniprot	P06683
GeneID	12279;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.172ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C9

Note: This product is for in vitro research use only