

Human C6/Complement component C6 ELISA Kit

Catalog No: #EK19682

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human C6/Complement component C6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Complement component C6
Target Name	Complement component C6
Accession No.	P13671
Uniprot	P13671
GeneID	729;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C6

Note: This product is for in vitro research use only