

Human BCL11A/B-cell lymphoma/leukemia 11A ELISA Kit



Catalog No: #EK19683

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BCL11A/B-cell lymphoma/leukemia 11A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human B-cell lymphoma/leukemia 11A
Target Name	B-cell lymphoma/leukemia 11A
Other Names	BCL-11A,B-cell CLL/lymphoma 11A,COUP-TF-interacting protein 1,Ecotropic viral integration site 9 protein homolog,EVI-9,Zinc finger protein 856,CTIP1,EVI9,KIAA1809,ZNF856
Accession No.	Q9H165
Uniprot	Q9H165
GeneID	53335;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.178ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BCL11A

Note: This product is for in vitro research use only