

Mouse Bcl11a/B-cell lymphoma/leukemia 11A ELISA Kit



Catalog No: #EK19684

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bcl11a/B-cell lymphoma/leukemia 11A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell lymphoma/leukemia 11A
Target Name	B-cell lymphoma/leukemia 11A
Other Names	BCL-11A,B-cell CLL/lymphoma 11A,COUP-TF-interacting protein 1,Ecotropic viral integration site 9 protein,EVI-9,Ctip1,Evi9,Kiaa1809
Accession No.	Q9QYE3
Uniprot	Q9QYE3
GeneID	14025;

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.312ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bcl11a

Note: This product is for in vitro research use only