

Mouse Ca3/Carbonic anhydrase 3 ELISA Kit

Catalog No: #EK19686

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ca3/Carbonic anhydrase 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carbonic anhydrase 3
Target Name	Carbonic anhydrase 3
Other Names	Carbonate dehydratase III, Carbonic anhydrase III, CA-III, Car3
Accession No.	P16015
Uniprot	P16015
GeneID	12350;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.061ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ca3

Note: This product is for in vitro research use only