

# Mouse Amfr/E3 ubiquitin-protein ligase AMFR ELISA Kit



Catalog No: #EK19687

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Amfr/E3 ubiquitin-protein ligase AMFR ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse E3 ubiquitin-protein ligase AMFR
Target Name	E3 ubiquitin-protein ligase AMFR
Other Names	Autocrine motility factor receptor,AMF receptor,RING-type E3 ubiquitin transferase AMFR
Accession No.	Q9R049
Uniprot	Q9R049
GeneID	23802;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.27ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Amfr

Note: This product is for in vitro research use only