

## Mouse Cd109/CD109 antigen ELISA Kit

Catalog No: #EK19688

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd109/CD109 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CD109 antigen
Target Name	CD109 antigen
Other Names	GPI-anchored alpha-2 macroglobulin-related protein,CD109
Accession No.	Q8R422
Uniprot	Q8R422
GeneID	235505;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name: Cd109

Note: This product is for in vitro research use only