

Rat Abca7/ATP-binding cassette sub-family A member 7 ELISA Kit



Catalog No: #EK19689

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Abca7/ATP-binding cassette sub-family A member 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ATP-binding cassette sub-family A member 7
Target Name	ATP-binding cassette sub-family A member 7
Accession No.	Q7TNJ2
Uniprot	Q7TNJ2
GeneID	299609;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Abca7

Note: This product is for in vitro research use only