

# Bovine ATP2A1/Sarcoplasmic/endoplasmic reticulum calcium ATPase 1 ELISA Kit



Catalog No: #EK19690

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ATP2A1/Sarcoplasmic/endoplasmic reticulum calcium ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Target Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Other Names	SERCA1, Calcium pump 1, Calcium-transporting ATPase sarcoplasmic reticulum type, fast twitch skeletal muscle isoform, Endoplasmic reticulum class 1/2 Ca(2+) ATPase
Accession No.	Q0VCY0
Uniprot	Q0VCY0
GeneID	518117;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ATP2A1

Note: This product is for in vitro research use only