

Rat Atp2a1/Sarcoplasmic/endoplasmic reticulum calcium ATPase 1 ELISA Kit

Catalog No: #EK19693

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Atp2a1/Sarcoplasmic/endoplasmic reticulum calcium ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Target Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Other Names	SERCA1, Calcium pump 1, Calcium-transporting ATPase sarcoplasmic reticulum type, fast twitch skeletal muscle isoform, Endoplasmic reticulum class 1/2 Ca(2+) ATPase
Accession No.	Q64578
Uniprot	Q64578
GeneID	116601;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.09ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp2a1

Note: This product is for in vitro research use only