

## Mouse Cntn1/Contactin-1 ELISA Kit

Catalog No: #EK19696

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cntn1/Contactin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Contactin-1
Target Name	Contactin-1
Other Names	Neural cell surface protein F3
Accession No.	P12960
Uniprot	P12960
GeneID	12805;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:31pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cntn1

Note: This product is for in vitro research use only