

Mouse Chad/Chondroadherin ELISA Kit

Catalog No: #EK19698

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Chad/Chondroadherin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Chondroadherin
Target Name	Chondroadherin
Other Names	Cartilage leucine-rich protein
Accession No.	O55226
Uniprot	O55226
GeneID	12643;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.162ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Chad

Note: This product is for in vitro research use only