

Rat Chad/Chondroadherin ELISA Kit

Catalog No: #EK19699

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Chad/Chondroadherin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chondroadherin
Target Name	Chondroadherin
Other Names	Cartilage leucine-rich protein
Accession No.	O70210
Uniprot	O70210
GeneID	29195;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.06ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Chad

Note: This product is for in vitro research use only