

Bovine COLEC11/Collectin-11 ELISA Kit

Catalog No: #EK19700



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bovine COLEC11/Collectin-11 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Collectin-11
Target Name	Collectin-11
Other Names	Collectin kidney protein 1,CL-K1
Accession No.	Q17QH6
Uniprot	Q17QH6
GeneID	613991;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.065ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:COLEC11

Note: This product is for in vitro research use only