

Mouse Colec11/Collectin-11 ELISA Kit

Catalog No: #EK19701

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Colec11/Collectin-11 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Collectin-11
Target Name	Collectin-11
Other Names	Collectin kidney protein 1,CL-K1
Accession No.	Q3SXB8
Uniprot	Q3SXB8
GeneID	71693;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.051ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Colec11

Note: This product is for in vitro research use only