

## Rat Atg5/Autophagy protein 5 ELISA Kit

Catalog No: #EK19702

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Atg5/Autophagy protein 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Autophagy protein 5
Target Name	Autophagy protein 5
Accession No.	Q3MQ06
Uniprot	Q3MQ06

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 32pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Atg5

Note: This product is for in vitro research use only