

Mouse Atp7a/Copper-transporting ATPase 1 ELISA Kit

Catalog No: #EK19703

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Atp7a/Copper-transporting ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Copper-transporting ATPase 1
Target Name	Copper-transporting ATPase 1
Other Names	Copper pump 1, Menkes disease-associated protein homolog, Mnk
Accession No.	Q64430
Uniprot	Q64430
GeneID	11977;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.092ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp7a

Note: This product is for in vitro research use only