

## Rat Atp7a/Copper-transporting ATPase 1 ELISA Kit

Catalog No: #EK19704

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Atp7a/Copper-transporting ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Copper-transporting ATPase 1
Target Name	Copper-transporting ATPase 1
Other Names	Copper pump 1, Menkes disease-associated protein homolog, Mnk
Accession No.	P70705
Uniprot	P70705
GeneID	24941;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.5pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Atp7a

Note: This product is for in vitro research use only