

# Human CPS1/Carbamoyl-phosphate synthase [ammonia], mitochondrial ELISA Kit



Catalog No: #EK19705

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CPS1/Carbamoyl-phosphate synthase [ammonia], mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carbamoyl-phosphate synthase [ammonia], mitochondrial
Target Name	Carbamoyl-phosphate synthase [ammonia], mitochondrial
Other Names	Carbamoyl-phosphate synthetase I,CPSase I
Accession No.	P31327
Uniprot	P31327
GeneID	1373;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:14.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CPS1

Note: This product is for in vitro research use only