

Mouse Cdk5r1/Cyclin-dependent kinase 5 activator 1 ELISA Kit



Catalog No: #EK19706

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Cdk5r1/Cyclin-dependent kinase 5 activator 1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cyclin-dependent kinase 5 activator 1 |
| Target Name | Cyclin-dependent kinase 5 activator 1 |
| Other Names | CDK5 activator 1, Cyclin-dependent kinase 5 regulatory subunit 1, TPKII regulatory subunit, Cdk5r, Nck5a |
| Accession No. | P61809 |
| Uniprot | P61809 |
| GeneID | 12569; |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdk5r1

Note: This product is for in vitro research use only