

Human ADTRP/Androgen-dependent TFPI-regulating protein ELISA Kit



Catalog No: #EK19708

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ADTRP/Androgen-dependent TFPI-regulating protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Androgen-dependent TFPI-regulating protein
Target Name	Androgen-dependent TFPI-regulating protein
Other Names	C6orf105
Accession No.	NP_116133.1/ncbi

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/ml Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADTRP

Note: This product is for in vitro research use only