

Bovine CASR/Extracellular calcium-sensing receptor ELISA Kit



Catalog No: #EK19710

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CASR/Extracellular calcium-sensing receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Extracellular calcium-sensing receptor
Target Name	Extracellular calcium-sensing receptor
Other Names	CaSR, Parathyroid cell calcium-sensing receptor, PCaR1, GPRC2A, PCAR1
Accession No.	P35384
Uniprot	P35384
GeneID	281038;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.176ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CASR

Note: This product is for in vitro research use only