

Mouse Casr/Extracellular calcium-sensing receptor ELISA Kit

Catalog No: #EK19712

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Casr/Extracellular calcium-sensing receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Extracellular calcium-sensing receptor
Target Name	Extracellular calcium-sensing receptor
Other Names	CaSR,Parathyroid cell calcium-sensing receptor,PCaR1,Gprc2a
Accession No.	Q9QY96
Uniprot	Q9QY96
GeneID	12374;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.09ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Casr

Note: This product is for in vitro research use only