

# Rat Casr/Extracellular calcium-sensing receptor ELISA Kit

Catalog No: #EK19713

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Casr/Extracellular calcium-sensing receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Extracellular calcium-sensing receptor
Target Name	Extracellular calcium-sensing receptor
Other Names	CaSR,Parathyroid cell calcium-sensing receptor,PCaR1,Gprc2a,Pcar1
Accession No.	P48442
Uniprot	P48442
GeneID	24247;

## Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.194ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Casr

Note: This product is for in vitro research use only