

# Human ATR/Serine/threonine-protein kinase ATR ELISA Kit



Catalog No: #EK19714

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ATR/Serine/threonine-protein kinase ATR ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Serine/threonine-protein kinase ATR
Target Name	Serine/threonine-protein kinase ATR
Other Names	Ataxia telangiectasia and Rad3-related protein,FRAP-related protein 1,FRP1
Accession No.	Q13535
Uniprot	Q13535
GeneID	545;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.066ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ATR

Note: This product is for in vitro research use only