

Human ASCL1/Achaete-scute homolog 1 ELISA Kit

Catalog No: #EK19715

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ASCL1/Achaete-scute homolog 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Achaete-scute homolog 1
Target Name	Achaete-scute homolog 1
Other Names	ASH-1, Class A basic helix-loop-helix protein 46, bHLHa46, ASH1, BHLHA46, HASH1
Accession No.	P50553
Uniprot	P50553
GeneID	429;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.66pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ASCL1

Note: This product is for in vitro research use only